Light path changeable in accordance with right or left better eye through one-touch operation

The light path to capture images into the digital camera can be changed in accordance with the examiner's better eye. By adjusting the dial on the back, you can change the light path through one-touch operation. The fundus can be observed without worrying about displacement between an image in the field of view and a captured image.



Yellow filter IN/OUT mechanism

High-contrast images can be obtained by using a built-in yellow filter in cornea fluorescence observation. Easy filter IN or OUT is available through one-touch operation. In its setting with our slitlamp, high-contrast images can be obtained by the best matching with a cobalt blue filter.

Camera models compatible with TD-2

TD-2 types do not come with a digital camera. For the purchase of a digital camera, refer to the following compatible cameras:

Compatible cameras manufactured by Canon EOS kiss Digital N

EOS DIGITAL REBEL XT/350D DIGITAL EOS DIGITAL REBEL XT/300D DIGITAL EOS 10D, EOS 20D, EOS D30, EOS D60



Capturing of still image with one push

By connecting a trigger cable between the slitlamp and camera, a still image can be captured every time the joystick button of the slitlamp is pushed.



Slitlamp models compatible with TD-2 types




Type A
Compatible models
SM-70N, SM-70,
SM-90N, SM-90



Type B Compatible models SM-2N, SM-2, SM-4N, SM-4

The digital camera adapter TD-2 allows you to import the images taken into the PC and display them on the monitor by installing a commercially available single-lens reflex digital camera and combining it with a slitlamp.





 \bullet Design and specifications are subject to change as improvements are made to the product.



