

Support precise and various kinds of test

System Networking

performance provides precise and

Fast switching speed of the lens

Convenient control board

Various convenience functions

A wide range of viewing angle



# PDR-7000

### SPECIFICATION

	measurement ra	-29.00 ~ +26.75D(general)	
	measurement	-19.00 ~ +16.75D(cross cylinder or prism)	
	measurement u	it 0.12D/0.25D/0.50D/1.00D	
	measurement ra	nge 0D ~ ±8.75D	
	measurement u	it 0.25D/1D/2D/3D	
	measurement ra	nge 0°~180°	
	measurement u	it 1*/5*/15*	
	measurement ra	nge 48 ~ 80mm, measurement distance, 35 ~ 70cm	
	measurement u	it 0.5mm/1mm	
	measurement ra	nge 0△~20△	
	measurement u	it 0.1△/0.5△/2△	
	measurement ra	nge 0°~360°	
	measurement u	it 1°/5°	
	automatic cross cylinder(±0.25D)		
	0.25D±cross cylinder		
	0.50D±cross cylinder		
	P.D/cover/pin hole(Ф2mm)/polarized light filter(45"/135")		
	madox road(right eye:horizontality, left eye:verticality)		
	red-green filter(right eye:red filter, left eye:green filter)		
	polarized light filter(right eye:135",45") (left eye:45",135")		
	seperation prizm(right eye:6△BU) (left eye:10△BU)		
	Refractor 323-389(W)x245(D)x59(H)mm/ 3.7kg		
	Controller 260(W)x230(D)x58(H)mm/ 3.1kg(printer included)		
	Junction box 210(W)x190(D)x70(H)mm/ 2.1kg		

SYSTEM NETWORKING



POTEC Co., Ltd.

1324, Gwanpyeong-dong, Yuseong-gu,
Daejeon 305-509, Kores
Tel +82 42 53 2336 Fax +82 42 632 3337
webmaster@potec.biz www.potec.biz

www. potec. biz

Make smart technology yours



PDR-7000



CE



# ontometry

Supply dual cross cylinder lens

Dual cross cylinder lens makes more accurate and easy astigmatism and proof test.

Automatic cover function

Since when reversing lens and switching mode, automatic cover function prevents intervention of accommodation force, easy and correct test is possible.

Congestion Contro

Test can be done by the center of the lens eye test in window both presbyopia and nearsighted test.

Nearsightedness test of all kinds

Various target makes possible many kinds of nearsightedness test.

Various test method

PDR-7000 provides various muscle balance test methods such as Von Graefe Test, Schober Test, Madox Rod Test, etc.

# design

Ergonomic design, provide wide viewing angle of 40-degree

Compact and thin size makes you test more conveniently because it covers small area of patient's face.

Test results can be secured practically and accurately due to wide range of viewing angle of 40 degree.

Forehead pedestal sensor

Forehead location detection, equipped with alert features such as indicating optimum face location for test, can check whether patients maintain correct posture or not.

Left / Right Independent PD Adjustment and fast speed of lens switching

As Left / Right PD adjusted independently, it enhances the reliability of test results and when using an eye chart, enough illumination can be obtained by using light. In addition, as the fast lens switching speed minimizes eye control and eye fatigue that can be migrated during a test, the quality of optometry is improved.

Provide easy maintenance and upgrade

As the forehead pedestal, facial contact, and eyepiece protect lens which are designed detachable type can be removed and cleaned easily, you can always maintain hygiene test environment. In addition, the software can be upgraded easily through USB.

# PDR-7000







Touch Screen Monitor

Clicking Jog Dial

Vision Test







# interface

#### Provide Tilting & Swivel monitor

Regardless of the position of the patient and the optometrist, the screen can be shared any angle and direction.

## Provide large screen and convenient touch screen display

As progress of diagnose and figures are expressed in an 10.4 inch wide screen with intuitive interface and various colors, it is easy to understand and aware of the test contents, and due to adopting touch-screen, optometry progress can be progressed easily even if contemptative due of families with the operations.

### Jog dial of click function and improved tact switch

Click function provides lens conversion and program execution much faster and conveniently and the voluminous keys arranged by function provide good touch sensitivity and improved durability.

#### Useful and convenient Vision Test features

Due to the useful Vision Test function, it is easy to explain to patients and make them understand about various refractive errors and diseases as well as near visual actiuity with near eye chart and auxiliary eye chart, and eye charts can be selected depending on user's taste and the purpose.



#### Drive automatically by one-touching the arm

#### 90 degree rotated table

POS-7000, providing optimum testing environment under any conditions, is a simple and convenient compact unit table.











# POTEC Optometry System

Potec's high-quality optometry system optimized for high-tech optometry equipment.

#### **SPECIFICATION**

1,240(W)x1,151(D)x1,875(H)mm
1,090(W)x480(D)x29(T)mm
90°
535mm(left-right)
±70mm
180mm
12V halogen
220V 5A
300W (Max)

M Design and specifications can be changed without prior notice for improveme

