


Air Injection & Irrigation Cannula

Disposable ophthalmic cannula used for maintaining or forming anterior chamber and irrigating cornea during surgery. It may be used for injection or removal of air, fluids & viscoelastic substances during an ophthalmic surgery.



Product Description

Product Code	Description with Gauge	Needle Length	Bend Length	Qty. per box
MIPL/E8/A5	Anterior Chamber Cannula (Rycroft) - 16G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A6	Anterior Chamber Cannula (Rycroft) - 19G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A7	Anterior Chamber Cannula (Rycroft) - 20G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A8	Anterior Chamber Cannula (Rycroft) - 21G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A1	Anterior Chamber Cannula (Rycroft) - 23G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A2	Anterior Chamber Cannula (Rycroft) - 25G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A9	Anterior Chamber Cannula (Rycroft) - 26G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A3	Anterior Chamber Cannula (Rycroft) - 27G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A10	Anterior Chamber Cannula (Rycroft) - 28G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A4	Anterior Chamber Cannula (Rycroft) - 30G	22 mm	6 mm	10 Pcs. (Sterile)

Product Code	Description with Gauge	Needle Length	Bend Length	Qty. per box
MIPL/E8/B1	Viscoelastic Cannula - 23G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E8/B2	Viscoelastic Cannula - 25G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E8/B3	Viscoelastic Cannula - 27G	22 mm	9 mm	10 Pcs. (Sterile)

Features

- Thin wall cannula used for maximum flow
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)