


Cystotome Cannula

Disposable ophthalmic cannula used for capsulorhexis, intercapsular/ endocapsular, and can opener capsulotomy technique.



Product Description

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E7/A1	Cystotome (Formed), 25G	16 mm	10 Pcs. (Sterile)
MIPL/E7/A5	Cystotome (Formed), 26G	16 mm	10 Pcs. (Sterile)
MIPL/E7/A2	Cystotome (Formed), 27G	16 mm	10 Pcs. (Sterile)
MIPL/E7/A3	Cystotome (Formed), 30G	16 mm	10 Pcs. (Sterile)
MIPL/E7/A4	Cystotome (Blumenthal), 27G	16 mm	10 Pcs. (Sterile)



Features

- Thin wall cannula used for maximum flow
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)