

PACHPEN

Equipped with



*Precision at Your
Fingertips*

Accurate - Achieve improved patient outcomes with Accutome's revolutionary digital technology.

IOP Correction Calculation - Converts IOP measurements in seconds.

Easy to Use - Superior ergonomic design fits comfortably in any hand.

Versatile - Helps to save time and ease your practice demands.

In 1982, we created the first ophthalmic pachymeter. Today, we are creating the best.



Accutome

263 Great Valley Parkway
Malvern, PA 19355 USA
800-979-2020 610-889-0200
fax 610-889-3233
www.accutome.com

THE NEXT GENERATION IN PACHYMETRY

**Gentle Touch Probe
ensures patient
comfort**

**Stores and averages
9 measurements**

**NO NEED TO
RECHARGE -
Long-lasting
Lithium Battery
will last approx.
15,000 readings**

Range of 300 to 999 microns
Accuracy of +/-5 microns
Resolution of +/- 1 micron
65 MHz sampling probe
Lightweight – 3 oz. (85 g.)


ACCUTOME®

Accutome
263 Great Valley Parkway
Malvern, PA 19355 USA
800-979-2020 610-889-0200
Fax 610-889-3233
www.accutome.com
eyeopener@accutome.com